


Product Datasheet

CD62P (SELP) antibody (orb1319854)

Description	Rabbit polyclonal antibody to CD62P (SELP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	SELP
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16109
NCBI	NM_003005
Dilution Range	WB 1:20000,ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...