


Product Datasheet

CTBS antibody (orb1319869)

Description	Rabbit polyclonal antibody to CTBS
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	DNA immunization. This antibody was made against a protein fragment from the C Terminus Region
Target	CTBS
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01459
NCBI	NM_004388
Dilution Range	WB 1:20000, IHC 1:150, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...