


Product Datasheet

NAP1L3 antibody (orb1319873)

Description	Rabbit polyclonal antibody to NAP1L3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	NAP1L3
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99457
NCBI	NM_004538
Dilution Range	WB 1:20000, IHC 1:150, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...