


## Product Datasheet

### XPC antibody (orb1319875)

|                            |                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to XPC                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                       |
| <b>Reactivity</b>          | Human                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                 |
| <b>Tested Applications</b> | ELISA                                                                                        |
| <b>Immunogen</b>           | DNA immunization. This antibody is specific for the N Terminus Region of the target protein. |
| <b>Target</b>              | XPC                                                                                          |
| <b>Note</b>                | For research use only                                                                        |
| <b>Isotype</b>             | IgG                                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q01831</a>                                                                       |
| <b>NCBI</b>                | <a href="#">NM_004628</a>                                                                    |
| <b>Dilution Range</b>      | ELISA: 1:100-1:2000                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                              |



ELISA using  
XPC  
antibody