


Product Datasheet

OLIG2 antibody (orb1319895)

Description	Rabbit polyclonal antibody to OLIG2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	DNA immunization. This antibody was made against a protein fragment from the N Terminus Region
Target	OLIG2
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13516
NCBI	NM_005806
Dilution Range	WB 1:20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...