


Product Datasheet

TCIRG1 antibody (orb1319900)

Description	Rabbit polyclonal antibody to TCIRG1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	TCIRG1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13488
NCBI	NM_006019
Dilution Range	WB 1:5000~20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...