


Product Datasheet

Nestin (NES) antibody (orb1319907)

Description	Rabbit polyclonal antibody to Nestin (NES)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the C Terminus Region of the target protein.
Target	NES
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P48681
NCBI	NM_006617
Dilution Range	ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
Nestin (NES)
antibody