


Product Datasheet

PHF16 (JADE3) antibody (orb1319932)

Description	Rabbit polyclonal antibody to PHF16 (JADE3)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA immunization. This antibody is specific for the C Terminus Region of the target protein.
Target	JADE3
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92613
NCBI	NM_014735
Dilution Range	WB 1:5000~20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...