


Product Datasheet

EFHC1 antibody (orb1319949)

Description	Rabbit polyclonal antibody to EFHC1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	DNA immunization. This antibody is specific for the C Terminus Region of the target protein.
Target	EFHC1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5JVL4
NCBI	NM_018100
Dilution Range	WB 1:20000,ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...