


Product Datasheet

KIAA1199 (CEMIP) antibody (orb1319956)

Description	Rabbit polyclonal antibody to KIAA1199 (CEMIP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	CEMIP
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WUJ3
NCBI	NM_018689
Dilution Range	WB: 1:5000-1:20000; ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using KIAA1199 (CEMIP) antibody