


Product Datasheet

PEPP2 (PLEKHA5) antibody (orb1319960)

Description	Rabbit polyclonal antibody to PEPP2 (PLEKHA5)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	PLEKHA5
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HAU0
NCBI	NM_019012
Dilution Range	ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
PEPP2
(PLEKHA5)
antibody