


Product Datasheet

TGF beta induced factor 2 (TGIF2) antibody (orb1319967)

Description	Rabbit polyclonal antibody to TGF beta induced factor 2 (TGIF2)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	TGIF2
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9GZN2
NCBI	NM_021809
Dilution Range	WB 1:20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...