


Product Datasheet

RAB11FIP1 antibody (orb1319978)

Description	Rabbit polyclonal antibody to RAB11FIP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	DNA immunization. This antibody was made against a protein fragment from the N Terminus Region
Target	RAB11FIP1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6WKZ4
NCBI	NM_025151
Dilution Range	WB 1:5000~20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...