


Product Datasheet

Netrin G2 (NTNG2) antibody (orb1319983)

Description	Rabbit polyclonal antibody to Netrin G2 (NTNG2)	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA	ELISA using Netrin G2 (NTNG2) antibody
Immunogen	DNA immunization. This antibody is specific for the C Terminus Region of the target protein.	
Target	NTNG2	
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
Uniprot ID	Q96CW9	
NCBI	NM_032536	
Dilution Range	WB: 1:5000-1:20000; ELISA: 1:100-1:2000	
Expiration Date	12 months from date of receipt.	