


Product Datasheet

FCP1 (CTDP1) antibody (orb1319990)

Description	Rabbit polyclonal antibody to FCP1 (CTDP1)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	CTDP1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5B0
NCBI	NM_048368
Dilution Range	WB: 1:5000-1:20000; ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
FCP1
(CTDP1)
antibody