


Product Datasheet

Kallikrein 7 (KLK7) antibody (orb1319998)

Description	Rabbit polyclonal antibody to Kallikrein 7 (KLK7)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	KLK7
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49862
NCBI	NM_139277
Dilution Range	WB: 1:500; ELISA: 1:100-1:2000; IHC: 1:600; IF: 1:150
Expiration Date	12 months from date of receipt.



ELISA using Kallikrein 7 (KLK7) antibody