


Product Datasheet

PALM2AKAP2 antibody (orb1320003)

Description	Rabbit polyclonal antibody to PALM2AKAP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	PALM2AKAP2
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2D5
NCBI	NM_147150
Dilution Range	ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
PALM2AKAP2
antibody