

Product Datasheet

PPP2R3A antibody (orb1320035)

Description	Rabbit polyclonal antibody to PPP2R3A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	PPP2R3A
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q06190
NCBI	NM_181897
Dilution Range	WB: 1:5000-1:20000; ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
PPP2R3A
antibody