


Product Datasheet

Fibronectin (FN1) antibody (orb1320056)

Description	Rabbit polyclonal antibody to Fibronectin (FN1)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	FN1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02751
NCBI	NM_212474
Dilution Range	ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using Fibronectin (FN1) antibody