


Product Datasheet

PKC zeta (PRKCZ) antibody (orb1320160)

Description	Rabbit polyclonal antibody to PKC zeta (PRKCZ)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	PRKCZ
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q05513
NCBI	NM_002744
Dilution Range	WB 1:20000,ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...