


Product Datasheet

PDPK1 antibody (orb1320162)

Description	Rabbit polyclonal antibody to PDPK1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	PDPK1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15530
NCBI	NM_002613
Dilution Range	WB 1:20000,ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using
PDPK1
antibody