

Product Datasheet

ATF 4 (ATF4) antibody (orb1320182)

Description	Rabbit polyclonal antibody to ATF 4 (ATF4)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Immunogen	DNA immunization. This antibody was made against a protein fragment from the Middle Region
Target	ATF4
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18848
NCBI	NM_001675
Dilution Range	WB: 1:100-1:2000; ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.