


Product Datasheet

Chromogranin A (CHGA) antibody (orb1320185)

Description	Rabbit polyclonal antibody to Chromogranin A (CHGA)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	CHGA
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10645
NCBI	NM_001275
Dilution Range	ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using Chromogranin A (CHGA) antibo...