


## Product Datasheet

### C12orf26 (METTL25) antibody (orb1320256)

|                            |                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to C12orf26 (METTL25)                                              |
| <b>Species/Host</b>        | Mouse                                                                                        |
| <b>Reactivity</b>          | Human                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                 |
| <b>Tested Applications</b> | IHC, WB                                                                                      |
| <b>Immunogen</b>           | Full length human recombinant protein of human C12orf26(NP_115606) produced in HEK293T cell. |
| <b>Target</b>              | METTL25                                                                                      |
| <b>Note</b>                | For research use only                                                                        |
| <b>Isotype</b>             | IgG1                                                                                         |
| <b>Clonality</b>           | Monoclonal                                                                                   |
| <b>MW</b>                  | 68 kDa                                                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q8N6Q8</a>                                                                       |
| <b>NCBI</b>                | <a href="#">NM_032230</a>                                                                    |
| <b>Dilution Range</b>      | WB 1:4000, IHC 1:150                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                              |



Western blot analysis of HEK293T cells w...