


Product Datasheet

Sterol carrier protein 2 (SCP2) antibody (orb1320288)

Description	Mouse monoclonal antibody to Sterol carrier protein 2 (SCP2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SCP2(NP_002970) produced in HEK293T cell
Target	SCP2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	58.8 kDa
Uniprot ID	P22307
NCBI	NM_002979
Dilution Range	WB 1:400~4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...