


Product Datasheet

SCRN1 antibody (orb1320291)

Description	Mouse monoclonal antibody to SCRN1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SCRN1(NP_055581) produced in HEK293T cell.
Target	SCRN1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.2 kDa
Uniprot ID	Q12765
NCBI	NM_014766
Dilution Range	WB 1:400~4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...