


Product Datasheet

Angiomotin (AMOT) antibody (orb1320359)

Description	Mouse monoclonal antibody to Angiomotin (AMOT)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human AMOT(NP_573572) produced in HEK293T cell.
Target	AMOT
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	72.4 kDa
Uniprot ID	Q4VCS5
NCBI	NM_133265
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...