

Product Datasheet

PD-L1 (CD274) antibody (orb1320371)

Description	Mouse monoclonal antibody to PD-L1 (CD274)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB, xMAP
Immunogen	Full length human recombinant protein of human CD274(NP_054862) produced in HEK293T cell.
Target	CD274
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31 kDa
Uniprot ID	Q9NZQ7
NCBI	NM_014143
Dilution Range	WB 1:200~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...