

Product Datasheet

UCHL5IP (HAUS7) antibody (orb1320373)

Description Mouse monoclonal antibody to UCHL5IP (HAUS7)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, WB

Immunogen Full length human recombinant protein of human HAUS7(NP_059988) produced in HEK293T cell.

Target HAUS7

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 40.6 kDa

Uniprot ID [Q99871](#)

NCBI [NM_017518](#)

Dilution Range WB: 1:400-1:2000, IF: 1:50-1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...