


Product Datasheet

CEP68 antibody (orb1320389)

Description	Mouse monoclonal antibody to CEP68
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CEP68(NP_055962) produced in HEK293T cell.
Target	CEP68
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	80.9 kDa
Uniprot ID	Q76N32
NCBI	NM_015147
Dilution Range	WB 1:400~4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...