


Product Datasheet

MFAP3 antibody (orb1320392)

Description	Mouse monoclonal antibody to MFAP3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human MFAP3(NP_005918) produced in HEK293T cell.
Target	MFAP3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	40 kDa
Uniprot ID	P55082
NCBI	NM_005927
Dilution Range	WB 1:1000 -1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...