


Product Datasheet

CYBC1 antibody (orb1320406)

Description	Mouse monoclonal antibody to CYBC1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human C17orf62(NP_001028218) produced in HEK293T cell.
Target	CYBC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	20.6 kDa
Uniprot ID	Q9BQA9
NCBI	NM_001033046
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...