


Product Datasheet

CHAC1 antibody (orb1320412)

Description	Mouse monoclonal antibody to CHAC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human CHAC1(NP_077016) produced in HEK293T cell.
Target	CHAC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	24.2 kDa
Uniprot ID	Q9BUX1
NCBI	NM_024111
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...