


Product Datasheet

TXNDC (TMX1) antibody (orb1320446)

Description	Mouse monoclonal antibody to TXNDC (TMX1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TMX1(NP_110382) produced in HEK293T cell.
Target	TMX1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.6 kDa
Uniprot ID	Q9H3N1
NCBI	NM_030755
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...