


Product Datasheet

HRASLS3 (PLA2G16) antibody (orb1320470)

Description	Mouse monoclonal antibody to HRASLS3 (PLA2G16)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PLA2G16(NP_009000) produced in HEK293T cell.
Target	PLAAT3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17.8 kDa
Uniprot ID	P53816
NCBI	NM_007069
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...