


Product Datasheet

UPRT antibody (orb1320501)

Description	Mouse monoclonal antibody to UPRT
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human UPRT(NP_659489) produced in HEK293T cell.
Target	UPRT
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	33.6 kDa
Uniprot ID	Q96BW1
NCBI	NM_145052
Dilution Range	WB 1:1000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...