


Product Datasheet

POLR3GL antibody (orb1320503)

Description	Mouse monoclonal antibody to POLR3GL
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human POLR3GL(NP_115681) produced in HEK293T cell.
Target	POLR3GL
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9BT43
NCBI	NM_032305
Dilution Range	WB 1:4000, IHC 1:150, IF: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...