


Product Datasheet

IL6R antibody (orb1320508)

Description	Mouse monoclonal antibody to IL6R
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human IL6R(NP_000556) produced in HEK293T cell.
Target	IL6R
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	49.8 kDa
Uniprot ID	P08887
NCBI	NM_000565
Dilution Range	WB 1:4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...