


Product Datasheet

IL1 alpha (IL1A) antibody (orb1320514)

Description	Mouse monoclonal antibody to IL1 alpha (IL1A)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human IL1A(NP_000566) produced in HEK293T cell.
Target	IL1A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	30.4 kDa
Uniprot ID	P01583
NCBI	NM_000575
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (10ug)...