


Product Datasheet

CD44 antibody (orb1320575)

Description	Mouse monoclonal antibody to CD44
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human CD44(NP_000601) produced in HEK293T cell.
Target	CD44
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	79.2 kDa
Uniprot ID	P16070
NCBI	NM_000610
Dilution Range	WB 1:4000~8000, IHC 1:150, IF 1:100, FLOW 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...