


Product Datasheet

EPCAM antibody (orb1320596)

Description	Mouse monoclonal antibody to EPCAM
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human EPCAM(NP_002345) produced in HEK293T cell.
Target	EPCAM
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	32.7 kDa
Uniprot ID	P16422
NCBI	NM_002354
Dilution Range	WB 1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...