


Product Datasheet

TUBB1 antibody (orb1320600)

Description	Mouse monoclonal antibody to TUBB1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TUBB1(NP_110400) produced in HEK293T cell.
Target	TUBB1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	50.1 kDa
Uniprot ID	Q9H4B7
NCBI	NM_030773
Dilution Range	WB 1:4000, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...