


Product Datasheet

ACTH (POMC) antibody (orb1320621)

Description	Mouse monoclonal antibody to ACTH (POMC)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human POMC(NP_001030333) produced in HEK293T cell.
Target	POMC
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	26.6 kDa
Uniprot ID	P01189
NCBI	NM_001035256
Dilution Range	WB 1:4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...