


Product Datasheet

TUBB2B antibody (orb1320627)

Description	Mouse monoclonal antibody to TUBB2B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TUBB2B(NP_821080) produced in HEK293T cell.
Target	TUBB2B
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	49.8 kDa
Uniprot ID	Q9BVA1
NCBI	NM_178012
Dilution Range	WB 1:4000, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...