


Product Datasheet

CD163 antibody (orb1320692)

Description	Mouse monoclonal antibody to CD163
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human CD163(NP_004235) produced in HEK293T cell.
Target	CD163
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	125.3 kDa
Uniprot ID	Q86VB7
NCBI	NM_004244
Dilution Range	WB 1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...