


Product Datasheet

CD33 antibody (orb1320747)

Description	Mouse monoclonal antibody to CD33
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human CD33(NP_001763) produced in HEK293T cell.
Target	CD33
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	38 kDa
Uniprot ID	P20138
NCBI	NM_001772
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...