


Product Datasheet

RARRES1 antibody (orb1320851)

Description	Mouse monoclonal antibody to RARRES1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human RARRES1(NP_996846) produced in HEK293T cell.
Target	RARRES1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	33.1 kDa
Uniprot ID	P49788
NCBI	NM_206963
Dilution Range	WB 1:4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...