


Product Datasheet

MelanA (MLANA) antibody (orb1320877)

Description	Mouse monoclonal antibody to MelanA (MLANA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human MLANA(NP_005502) produced in HEK293T cell.
Target	MLANA
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	13 kDa
Uniprot ID	Q16655
NCBI	NM_005511
Dilution Range	WB 1:4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...