

Product Datasheet

Synaptophysin (SYP) antibody (orb1320881)

Description Mouse monoclonal antibody to Synaptophysin (SYP)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human SYP(NP_003170) produced in HEK293T cell.

Target SYP

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 33.7 kDa

Uniprot ID [P08247](#)

NCBI [NM_003179](#)

Dilution Range WB 1:4000, IHC 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...