


Product Datasheet

Parvin alpha (PARVA) antibody (orb1320917)

Description	Mouse monoclonal antibody to Parvin alpha (PARVA)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human PARVA(NP_060692) produced in HEK293T cell.
Target	PARVA
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	42.1 kDa
Uniprot ID	Q9NVD7
NCBI	NM_018222
Dilution Range	WB 1:200~4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...